

Anti-GREB1 antibody (STJA0004458)

STJA0004458

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-Protein Greb1 is suitable for use in Immunohistochemistry and ELISA research applications. |
| Applications | IHC/ELISA |
| Host/Source | Rabbit |
| Reactivity | Human/Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | Affinity purification |
| Dilution Range | IHC 1:50-1:200 |
| Formulation | PBS with 0.05% sodium azide and 50% glycerol, PH7.4 |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term, and avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|----------------------------------|
| Gene ID | 9687 |
| Gene Symbol | GREB1 |
| Uniprot ID | GREB1_HUMAN |
| Immunogen | Synthetic peptide of human GREB1 |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |