

Anti-ARHGAP15 antibody (STJA0004462)

STJA0004462

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Rho Gtpase-Activating Protein 15 is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
Applications	WB/IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.4 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:200-1:1000, IHC 1:50-1:200
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55843
Gene Symbol	ARHGAP15
Uniprot ID	RHG15_HUMAN
Immunogen	Recombinant protein of human ARHGAP15
Immunogen Region	
Specificity	
Immunogen Sequence	