

Anti-C9orf116 antibody (STJA0004464)

STJA0004464

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Piercer Of Microtubule Wall 1 Protein is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.96 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:50-1:300
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	138162
Gene Symbol	PIERCE1
Uniprot ID	PIRC1_HUMAN
Immunogen	Fusion protein of human C9orf116
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	