

Anti-ACSBG2 antibody (STJA0004465)

STJA0004465

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Long-Chain-Fatty-Acid--CoA Ligase Acsbg2 is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.3 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:100
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	81616
Gene Symbol	ACSBG2
Uniprot ID	ACBG2_HUMAN
Immunogen	Recombinant protein of human ACSBG2
Immunogen Region	
Specificity	
Immunogen Sequence	