

Anti-EEF1AKMT2 antibody (STJA0004471)

STJA0004471

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Eef1a Lysine Methyltransferase 2 is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.96 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:50-1:300
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	399818
Gene Symbol	EEF1AKMT2
Uniprot ID	EFMT2_HUMAN
Immunogen	Fusion protein of human EEF1AKMT2
Immunogen Region	
Specificity	
Immunogen Sequence	