

Anti-CCDC181 antibody (STJA0004486)

STJA0004486

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Coiled-Coil Domain-Containing Protein 181 is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.8 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:40-1:200
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57821
Gene Symbol	CCDC181
Uniprot ID	CC181_HUMAN
Immunogen	Fusion protein of human CCDC181
Immunogen Region	
Specificity	
Immunogen Sequence	