

Anti-GKAP1 antibody (STJA0004488)

STJA0004488

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-G Kinase-Anchoring Protein 1 is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.62 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:90-1:450
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	80318
Gene Symbol	GKAP1
Uniprot ID	GKAP1_HUMAN
Immunogen	Fusion protein of human GKAP1
Immunogen Region	
Specificity	
Immunogen Sequence	