

Anti-MEDAG antibody (STJA0004494)

STJA0004494

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-Mesenteric Estrogen-Dependent Adipogenesis Protein is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications. |
| Applications | WB/IHC/ELISA |
| Host/Source | Rabbit |
| Reactivity | Human/Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 1.2 mg/mL |
| Conjugation | Unconjugated |
| Purification | Antigen affinity purification |
| Dilution Range | WB 1:500-1:2000, IHC 1:50-1:300 |
| Formulation | PBS with 0.05% Na ₃ N and 40% Glycerol, pH7.4 |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term, and avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|----------------------------------|
| Gene ID | 84935 |
| Gene Symbol | MEDAG |
| Uniprot ID | MEDAG_HUMAN |
| Immunogen | Synthetic peptide of human MEDAG |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |