

Anti-CARD17 antibody (STJA0004496)

STJA0004496

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Putative Caspase Recruitment Domain-Containing Protein 17p is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.7 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:200
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	CARD17P
Uniprot ID	CAR17_HUMAN
Immunogen	Synthetic peptide of human CARD17
Immunogen Region	
Specificity	
Immunogen Sequence	