

**Anti-CYB5D1 antibody (STJA0004508)**

STJA0004508

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Cytochrome B5 Domain-Containing Protein 1 is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
<b>Applications</b>	WB/IHC/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1.4 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Antigen affinity purification
<b>Dilution Range</b>	WB 1:500-1:2000, IHC 1:40-1:200
<b>Formulation</b>	PBS with 0.05% NaN <sub>3</sub> and 40% Glycerol, pH7.4
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term, and avoid freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	124637
<b>Gene Symbol</b>	CYB5D1
<b>Uniprot ID</b>	CB5D1_HUMAN
<b>Immunogen</b>	Synthetic peptide of human CYB5D1
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	