

## Anti-GLIPR1L1 antibody (STJA0004515)

STJA0004515

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Glipr1-Like Protein 1 is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
<b>Applications</b>	WB/IHC/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.6 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Antigen affinity purification
<b>Dilution Range</b>	WB 1:500-1:2000, IHC 1:25-1:100
<b>Formulation</b>	PBS with 0.05% Na <sub>3</sub> N and 40% Glycerol, pH7.4
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term, and avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	256710
<b>Gene Symbol</b>	GLIPR1L1
<b>Uniprot ID</b>	GPRL1_HUMAN
<b>Immunogen</b>	Fusion protein of human GLIPR1L1
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	