

Anti-HSDL2 antibody (STJA0004527)

STJA0004527

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Hydroxysteroid Dehydrogenase-Like Protein 2 is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.38 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:50-1:200
Formulation	PBS with 0.05% NaN3 and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84263
Gene Symbol	HSDL2
Uniprot ID	HSDL2_HUMAN
Immunogen	Fusion protein of human HSDL2
Immunogen Region	
Specificity	
Immunogen Sequence	