

Anti-OTUD6A antibody (STJA0004542)

STJA0004542

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Otu Domain-Containing Protein 6a is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.44 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	139562
Gene Symbol	OTUD6A
Uniprot ID	OTU6A_HUMAN
Immunogen	Synthetic peptide of human OTUD6A
Immunogen Region	
Specificity	
Immunogen Sequence	