

## Anti-COMMD6 antibody (STJA0004547)

STJA0004547

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Comm Domain-Containing Protein 6 is suitable for use in Immunohistochemistry and ELISA research applications.
<b>Applications</b>	IHC/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1.4 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Antigen affinity purification
<b>Dilution Range</b>	IHC 1:30-1:150
<b>Formulation</b>	PBS with 0.05% NaN3 and 40% Glycerol, pH7.4
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term, and avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">170622</a>
<b>Gene Symbol</b>	<a href="#">COMMD6</a>
<b>Uniprot ID</b>	<a href="#">COMMD6_HUMAN</a>
<b>Immunogen</b>	Full length fusion protein
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081