

Anti-TMED4 antibody (STJA0004551)

STJA0004551

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-Transmembrane Emp24 Domain-Containing Protein 4 is suitable for use in Immunohistochemistry and ELISA research applications. |
| Applications | IHC/ELISA |
| Host/Source | Rabbit |
| Reactivity | Human/Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.3 mg/mL |
| Conjugation | Unconjugated |
| Purification | Affinity purification |
| Dilution Range | IHC 1:30-150 |
| Formulation | PBS with 0.05% sodium azide and 50% glycerol, PH7.4 |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term, and avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|-------------------------------|
| Gene ID | 222068 |
| Gene Symbol | TMED4 |
| Uniprot ID | TMED4_HUMAN |
| Immunogen | Fusion protein of human TMED4 |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |