

**Anti-PRUNE antibody (STJA0004552)**

STJA0004552

**GENERAL INFORMATION**

|                          |                                                                                                                                        |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                                     |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-Exopolyphosphatase Prune1 is suitable for use in Immunohistochemistry and ELISA research applications. |
| <b>Applications</b>      | IHC/ELISA                                                                                                                              |
| <b>Host/Source</b>       | Rabbit                                                                                                                                 |
| <b>Reactivity</b>        | Human                                                                                                                                  |

**PRODUCT PROPERTIES**

|                            |                                                             |
|----------------------------|-------------------------------------------------------------|
| <b>Clonality</b>           | Polyclonal                                                  |
| <b>Clone ID</b>            |                                                             |
| <b>Concentration</b>       | 0.6 mg/mL                                                   |
| <b>Conjugation</b>         | Unconjugated                                                |
| <b>Purification</b>        | Antigen affinity purification                               |
| <b>Dilution Range</b>      | IHC 1:30-1:150                                              |
| <b>Formulation</b>         | PBS with 0.05% NaN <sub>3</sub> and 40% Glycerol, pH7.4     |
| <b>Isotype</b>             | IgG                                                         |
| <b>Storage Instruction</b> | Store at -20°C for long term, and avoid freeze-thaw cycles. |

**TARGET INFORMATION**

|                           |                            |
|---------------------------|----------------------------|
| <b>Gene ID</b>            | 58497                      |
| <b>Gene Symbol</b>        | PRUNE1                     |
| <b>Uniprot ID</b>         | PRUN1_HUMAN                |
| <b>Immunogen</b>          | Full length fusion protein |
| <b>Immunogen Region</b>   |                            |
| <b>Specificity</b>        |                            |
| <b>Immunogen Sequence</b> |                            |