

Anti-FOXRED2 antibody (STJA0004572)

STJA0004572

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Fad-Dependent Oxidoreductase Domain-Containing Protein 2 is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.66 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:40-1:200
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	80020
Gene Symbol	FOXRED2
Uniprot ID	FXRD2_HUMAN
Immunogen	Fusion protein of human FOXRED2
Immunogen Region	
Immunogen Specificity	
Immunogen Sequence	