

Anti-ADHFE1 antibody (STJA0004577)

STJA0004577

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Hydroxyacid-Oxoacid Transhydrogenase-Mitochondrial is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.3 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:200
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	137872
Gene Symbol	ADHFE1
Uniprot ID	HOT_HUMAN
Immunogen	Synthetic peptide of human ADHFE1
Immunogen Region	
Specificity	
Immunogen Sequence	