

Anti-UNC45B antibody (STJA0004578)

STJA0004578

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Protein Unc-45 Homolog B is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.8 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	146862
Gene Symbol	UNC45B
Uniprot ID	UN45B_HUMAN
Immunogen	Fusion protein of human UNC45B
Immunogen Region	
Specificity	
Immunogen Sequence	