

Anti-GSPT2 antibody (STJA0004583)

STJA0004583

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Eukaryotic Peptide Chain Release Factor Gtp-Binding Subunit Erf3b is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
Applications	WB/IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.78 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	WB 1:500-1:2000, IHC 1:40-1:200
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23708
Gene Symbol	GSPT2
Uniprot ID	ERF3B_HUMAN
Immunogen	Fusion protein of human GSPT2
Immunogen Region	
Specificity	
Immunogen Sequence	