

Anti-DYNAP antibody (STJA0004587)

STJA0004587

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Dynactin-Associated Protein is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.6 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:30-1:150
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	284254
Gene Symbol	DYNAP
Uniprot ID	DYNAP_HUMAN
Immunogen	Synthetic peptide of human DYNAP
Immunogen Region	
Specificity	
Immunogen	
Sequence	