

Anti-DHRSX antibody (STJA0004596)

STJA0004596

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Dehydrogenase/Reductase Sdr Family Member On Chromosome X is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
Applications	WB/IHC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.9 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	WB 1:500-1:2000, IHC 1:25-1:100
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	207063
Gene Symbol	DHRSX
Uniprot ID	DHRSX_HUMAN
Immunogen	Fusion protein of human DHRSX
Immunogen Region	
Specificity	
Immunogen Sequence	