

Anti-NBEA antibody (STJA0004614)

STJA0004614

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Neurobeachin is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.4 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:25-1:100
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	26960
Gene Symbol	NBEA
Uniprot ID	NBEA_HUMAN
Immunogen	Synthetic peptide of human NBEA
Immunogen Region	
Specificity	
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081