

Anti-HMGCLL1 antibody (STJA0004619)

STJA0004619

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-3-Hydroxy-3-Methylglutaryl-CoA Lyase-Cytoplasmic is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.02 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:50-1:300
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54511
Gene Symbol	HMGCLL1
Uniprot ID	HMGCLL1_HUMAN
Immunogen	Fusion protein of human HMGCLL1
Immunogen Region	
Specificity	
Immunogen Sequence	