

Anti-RAB3IL1 antibody (STJA0004622)

STJA0004622

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Guanine Nucleotide Exchange Factor For Rab-3a is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.32 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	WB 1:1000-1:5000
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5866
Gene Symbol	RAB3IL1
Uniprot ID	R3GEF_HUMAN
Immunogen	Fusion protein of human RAB3IL1
Immunogen Region	
Specificity	
Immunogen Sequence	