

Anti-SYTL5 antibody (STJA0004625)

STJA0004625

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Synaptotagmin-Like Protein 5 is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.7 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:25-1:100
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	94122
Gene Symbol	SYTL5
Uniprot ID	SYTL5_HUMAN
Immunogen	Recombinant protein of human SYTL5
Immunogen Region	
Specificity	
Immunogen Sequence	