

Anti-TBCK antibody (STJA0004627)

STJA0004627

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Tbc Domain-Containing Protein Kinase-Like Protein is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.78 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:30-1:150
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	93627
Gene Symbol	TBCK
Uniprot ID	TBCK_HUMAN
Immunogen	Synthetic peptide of human TBCK
Immunogen Region	
Specificity	
Immunogen Sequence	