

**Anti-TBCK antibody (STJA0004628)**

STJA0004628

**GENERAL INFORMATION**

|                          |                                                                                                                                                                |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                                                             |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-Tbc Domain-Containing Protein Kinase-Like Protein is suitable for use in Immunohistochemistry and ELISA research applications. |
| <b>Applications</b>      | IHC/ELISA                                                                                                                                                      |
| <b>Host/Source</b>       | Rabbit                                                                                                                                                         |
| <b>Reactivity</b>        | Human/Mouse                                                                                                                                                    |

**PRODUCT PROPERTIES**

|                            |                                                             |
|----------------------------|-------------------------------------------------------------|
| <b>Clonality</b>           | Polyclonal                                                  |
| <b>Clone ID</b>            |                                                             |
| <b>Concentration</b>       | 1.14 mg/mL                                                  |
| <b>Conjugation</b>         | Unconjugated                                                |
| <b>Purification</b>        | Antigen affinity purification                               |
| <b>Dilution Range</b>      | IHC 1:50-1:300                                              |
| <b>Formulation</b>         | PBS with 0.05% NaN <sub>3</sub> and 40% Glycerol, pH7.4     |
| <b>Isotype</b>             | IgG                                                         |
| <b>Storage Instruction</b> | Store at -20°C for long term, and avoid freeze-thaw cycles. |

**TARGET INFORMATION**

|                           |                                 |
|---------------------------|---------------------------------|
| <b>Gene ID</b>            | 93627                           |
| <b>Gene Symbol</b>        | TBCK                            |
| <b>Uniprot ID</b>         | TBCK_HUMAN                      |
| <b>Immunogen</b>          | Synthetic peptide of human TBCK |
| <b>Immunogen Region</b>   |                                 |
| <b>Specificity</b>        |                                 |
| <b>Immunogen Sequence</b> |                                 |