

Anti-TSPEAR antibody (STJA0004631)

STJA0004631

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Thrombospondin-Type Laminin G Domain And Ear Repeat-Containing Protein is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.3 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:200
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54084
Gene Symbol	TSPEAR
Uniprot ID	TSEAR_HUMAN
Immunogen	Recombinant protein of human TSPEAR
Immunogen Region	
Specificity	
Immunogen Sequence	