

**Anti-HMGB4 antibody (STJA0004636)**

STJA0004636

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-High Mobility Group Protein B4 is suitable for use in Western Blot and ELISA research applications.
<b>Applications</b>	WB/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.2 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB 1:200-1:1000
<b>Formulation</b>	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term, and avoid freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	127540
<b>Gene Symbol</b>	HMGB4
<b>Uniprot ID</b>	HMGB4_HUMAN
<b>Immunogen</b>	Recombinant protein of human HMGB4
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	