

## Anti-PPCDC antibody (STJA0004650)

STJA0004650

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Phosphopantothencycysteine Decarboxylase is suitable for use in Immunohistochemistry and ELISA research applications.
<b>Applications</b>	IHC/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.96 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Antigen affinity purification
<b>Dilution Range</b>	IHC 1:50-1:300
<b>Formulation</b>	PBS with 0.05% NaN <sub>3</sub> and 40% Glycerol, pH7.4
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term, and avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">60490</a>
<b>Gene Symbol</b>	<a href="#">PPCDC</a>
<b>Uniprot ID</b>	<a href="#">COAC_HUMAN</a>
<b>Immunogen</b>	Full length fusion protein
<b>Immunogen Region</b>	
<b>Immunogen Specificity</b>	
<b>Immunogen Sequence</b>	