

Anti-ACSF2 antibody (STJA0004652)

STJA0004652

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Medium-Chain Acyl-CoA Ligase Acsf2-Mitochondrial is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
Applications	WB/IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.96 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	WB 1:500-1:2000, IHC 1:30-1:150
Formulation	PBS with 0.05% NaN3 and 40% Glycerol, pH7.4
Isotype	IgG
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	80221
Gene Symbol	ACSF2
Uniprot ID	ACSF2_HUMAN
Immunogen	Fusion protein of human ACSF2
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	