

Anti-RAB3C antibody (STJA0004655)

STJA0004655

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-Ras-Related Protein Rab-3c is suitable for use in Western Blot and ELISA research applications. |
| Applications | WB/ELISA |
| Host/Source | Rabbit |
| Reactivity | Human/Mouse/Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.3 mg/mL |
| Conjugation | Unconjugated |
| Purification | Affinity purification |
| Dilution Range | WB 1:500-1:2000 |
| Formulation | PBS with 0.05% sodium azide and 50% glycerol, PH7.4 |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term, and avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|----------------------------------|
| Gene ID | 115827 |
| Gene Symbol | RAB3C |
| Uniprot ID | RAB3C_HUMAN |
| Immunogen | Synthetic peptide of human RAB3c |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |