

Anti-FAM84A antibody (STJA0004679)

STJA0004679

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Protein Lratd1 is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
Applications	WB/IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.02 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	WB 1:500-1:2000, IHC 1:40-1:200
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	151354
Gene Symbol	LRATD1
Uniprot ID	LRAT1_HUMAN
Immunogen	Synthetic peptide of human FAM84A
Immunogen Region	
Specificity	
Immunogen	
Sequence	