

Anti-AGAP3 antibody (STJA0004692)

STJA0004692

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Arf-Gap With Gtpase-Ank Repeat And Ph Domain-Containing Protein 3 is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.72 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:50-1:100
Formulation	PBS with 0.05% Na ₃ N and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	116988
Gene Symbol	AGAP3
Uniprot ID	AGAP3_HUMAN
Immunogen	Synthetic peptide of human AGAP3
Immunogen Region	
Immunogen Specificity	
Immunogen Sequence	