

Anti-KCNA7 antibody (STJA0004702)

STJA0004702

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Potassium Voltage-Gated Channel Subfamily A Member 7 is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
Applications	WB/IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.6 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000, IHC 1:50-1:200
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
Isotype	IgG
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3743
Gene Symbol	KCNA7
Uniprot ID	KCNA7_HUMAN
Immunogen	Synthetic peptide of human KCNA7
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	