

Anti-CCT8L2 antibody (STJA0004705)

STJA0004705

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-T-Complex Protein 1 Subunit Theta-Like 2 is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.7 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:30-1:150
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	150160
Gene Symbol	CCT8L2
Uniprot ID	TCPQM_HUMAN
Immunogen	Fusion protein of human CCT8L2
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	