

## Anti-FAM213A antibody (STJA0004722)

STJA0004722

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Peroxiredoxin-Like 2a is suitable for use in Western Blot and ELISA research applications.
<b>Applications</b>	WB/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse/Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1.02 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Antigen affinity purification
<b>Dilution Range</b>	WB 1:500-1:2000
<b>Formulation</b>	PBS with 0.05% NaN <sub>3</sub> and 40% Glycerol, pH7.4
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term, and avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">84293</a>
<b>Gene Symbol</b>	<a href="#">PRXL2A</a>
<b>Uniprot ID</b>	<a href="#">PXL2A_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide of human FAM213A
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	