

Anti-CCDC106 antibody (STJA0004734)

STJA0004734

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-Coiled-Coil Domain-Containing Protein 106 is suitable for use in Immunohistochemistry and ELISA research applications. |
| Applications | IHC/ELISA |
| Host/Source | Rabbit |
| Reactivity | Human/Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.5 mg/mL |
| Conjugation | Unconjugated |
| Purification | Affinity purification |
| Dilution Range | IHC 1:50-1:200 |
| Formulation | PBS with 0.05% sodium azide and 50% glycerol, PH7.4 |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term, and avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--------------------------------------|
| Gene ID | 29903 |
| Gene Symbol | CCDC106 |
| Uniprot ID | CC106_HUMAN |
| Immunogen | Recombinant protein of human CCDC106 |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |