

Anti-CDADC1 antibody (STJA0004738)

STJA0004738

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cytidine And Dcnp Deaminase Domain-Containing Protein 1 is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.4 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:30-1:150
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	81602
Gene Symbol	CDADC1
Uniprot ID	CDAC1_HUMAN
Immunogen	Full length fusion protein
Immunogen Region	
Specificity	
Immunogen Sequence	