

Anti-QTRT1 antibody (STJA0004742)

STJA0004742

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Queuine Trna-Ribosyltransferase Catalytic Subunit 1 is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.08 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	81890
Gene Symbol	QTRT1
Uniprot ID	TGT_HUMAN
Immunogen	Fusion protein of human QTRT1
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	