

Anti-CECR5 antibody (STJA0004743)

STJA0004743

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Haloacid Dehalogenase-Like Hydrolase Domain-Containing 5 is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.3 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27440
Gene Symbol	HDHD5
Uniprot ID	HDHD5_HUMAN
Immunogen	Synthetic peptide of human CECR5
Immunogen Region	
Specificity	
Immunogen Sequence	