

Anti-TRIM34 antibody (STJA0004746)

STJA0004746

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Tripartite Motif-Containing Protein 34 is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
Applications	WB/IHC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.4 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:200-1:1000, IHC 1:25-1:100
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	53840
Gene Symbol	TRIM34
Uniprot ID	TRI34_HUMAN
Immunogen	Synthetic peptide of human TRIM34
Immunogen Region	
Specificity	
Immunogen Sequence	