

**Anti-EPT1 antibody (STJA0004749)**

STJA0004749

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Ethanolaminephosphotransferase 1 is suitable for use in Western Blot research applications.
<b>Applications</b>	WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Antigen Affinity Purification
<b>Dilution Range</b>	WB: 1:1000-2000; IP: 3ug/sample; IF: 1:100-200
<b>Formulation</b>	PBS with 0.1% sodium azide, 50% glycerol, pH 7.4
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term, and avoid freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">85465</a>
<b>Gene Symbol</b>	<a href="#">SELENOI</a>
<b>Uniprot ID</b>	<a href="#">EPT1_HUMAN</a>
<b>Immunogen</b>	Recombinant human EPT1 protein expressed by E.coli
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081