

Anti-CARD9 antibody (STJA0004760)

STJA0004760

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Caspase Recruitment Domain-Containing Protein 9 is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.8 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:100-1:300
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64170
Gene Symbol	CARD9
Uniprot ID	CARD9_HUMAN
Immunogen	Recombinant protein of human CARD9
Immunogen Region	
Specificity	
Immunogen Sequence	