

## Anti-SMOC2 antibody (STJA0004766)

STJA0004766

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Sparc-Related Modular Calcium-Binding Protein 2 is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
<b>Applications</b>	WB/IHC/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB 1:200-1:1000, IHC 1:25-1:100
<b>Formulation</b>	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term, and avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	64094
<b>Gene Symbol</b>	SMOC2
<b>Uniprot ID</b>	SMOC2_HUMAN
<b>Immunogen</b>	Synthetic peptide of human SMOC2
<b>Immunogen Region</b>	
<b>Immunogen Specificity</b>	
<b>Immunogen Sequence</b>	