

Anti-PIEZO2 antibody (STJA0004767)

STJA0004767

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Piezo-Type Mechanosensitive Ion Channel Component 2 is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.62 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:50-1:300
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	63895
Gene Symbol	PIEZO2
Uniprot ID	PIEZ2_HUMAN
Immunogen	Synthetic peptide of human PIEZO2
Immunogen Region	
Specificity	
Immunogen Sequence	