

**Anti-ELAC1 antibody (STJA0004779)**

STJA0004779

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Zinc Phosphodiesterase Elac Protein 1 is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
<b>Applications</b>	WB/IHC/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.9 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Antigen affinity purification
<b>Dilution Range</b>	WB 1:1000-1:5000, IHC 1:50-1:300
<b>Formulation</b>	PBS with 0.05% NaN3 and 40% Glycerol, pH7.4
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term, and avoid freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">55520</a>
<b>Gene Symbol</b>	<a href="#">ELAC1</a>
<b>Uniprot ID</b>	<a href="#">RNZ1_HUMAN</a>
<b>Immunogen</b>	Fusion protein of human ELAC1
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	