

Anti-METTL17 antibody (STJA0004781)

STJA0004781

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Methyltransferase-Like Protein 17-Mitochondrial is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
Applications	WB/IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.5 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	WB 1:500-1:2000, IHC 1:20-1:100
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	64745
Gene Symbol	METTL17
Uniprot ID	MET17_HUMAN
Immunogen	Synthetic peptide of human METTL17
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	