

Anti-PLEKHF2 antibody (STJA0004785)

STJA0004785

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Pleckstrin Homology Domain-Containing Family F Member 2 is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
Applications	WB/IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.62 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	WB 1:1000-1:5000, IHC 1:100-1:300
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79666
Gene Symbol	PLEKHF2
Uniprot ID	PKHF2_HUMAN
Immunogen	Fusion protein of human PLEKHF2
Immunogen Region	
Specificity	
Immunogen Sequence	